

INSTRUCTIONS FOR USE

1. Check the condition of the ultrasonic cleaner's power plug and cord.
2. Pour the cleaning solution into the stainless steel tank. Make sure the fluid level remains between 1/3 and 2/3 of the height of the tank.
3. Connect the device to the mains and turn on the switch (I).
4. For appliances with heating function, select the temperature desired temperature by pressing the \uparrow or \downarrow arrows on the control panel. To start the heating cycle, press the «Heating» button (heating). A flashing dot on the temperature display indicates that the heating function is active. The liquid temperature is displayed in real time and adjusted to reach the desired temperature. Due to the selected temperature, the heating time may be long before the setpoint is reached.
5. Place the tools to be cleaned in the basket and then place the basket in the tub, leaning on the edges.
6. Select the running time by pressing the \uparrow or \downarrow arrows on the control panel. control (pressing the average between 3 and 30 minutes).
7. Press the «Ultrasound» button. The flashing dot on the minute display, the Ripples on the surface of the liquid and the sound of vibrations indicate that the cleaning cycle is in progress.
8. An acoustic signal indicates the end of the cleaning cycle.
9. Remove the basket from the tank and rinse the tools.
10. Empty the tank manually by turning the drain valve.

Precautions for use

1. Read the instructions carefully before use.
2. Check the electrical connection before use.
3. Keep the appliance away from any flame.
4. Keep flammable liquids away from the appliance.
5. Do not activate the ultrasonic or heating function if the tank does not contain cleaning liquid, so as not to damage the device.
6. Leave the filter firmly in place at the bottom of the tank.

Maintenance

1. Periodically check the general condition of the device.
2. The device must be used by qualified personnel.
3. Check the condition of the electrical connections, the tank and the drain valve.
4. Do not let dust accumulate on the fan grille of the ultrasonic emitters.

Contents

- 1 Ultrasonic bath
- 1 Lid
- 1 Power cord
- 1 Basket
- 1 User Manual



USER MANUAL



Technical characteristics

Model	Capacity	Frequency (KHz)	Ultrasound power (IN)	Heating power (IN)	Timer (Min)	Temperature (°C)
MHC 30U	3L	40	120	150	1-99	20 to 80
MHC 50U	5L			150		
MHC 100U	10L		240	350		
MHC 130U	13L		300	450		
MHC 190U	19L			450		



INSTRUCTION FOR USE

1. Check the condition of the plug and the power cable of the ultrasonic bath.
2. Pour the cleaning solution into the stainless steel tank. Make sure the fluid level remains between 1/3 and 2/3 the height of the tank.
3. Connect the device to the mains and turn on the switch (I).
4. For appliances equipped with the heating function, select the desired temperature by pressing the \uparrow or \downarrow arrows on the control panel. To start the heating cycle, press the «Heating» button (heating). A flashing point on the temperature display indicates that the heating function is on. The liquid temperature is displayed in real time and regulated to achieve the desired temperature. Following temperature selected, the heating time may be long before reaching the set point.
5. Place the instruments to be cleaned in the basket then place the basket in the tank, resting on the edges.
6. Select the operating time by pressing the \uparrow or \downarrow arrows on the control panel (by pressing the average between 3 and 30min).
7. Press the “ultrasound” button. The flashing dot on the minute display, the ripples on the surface of the liquid and the sound of vibrations indicate that the cleaning cycle is in progress.
8. An audible signal indicates the end of the cleaning cycle.
9. Remove the basket from the tank and rinse the instruments.
10. Empty the tank manually by turning the drain valve.

Precautions for use

1. Read the instructions carefully before use.
2. Check the electrical connection before use.
3. Keep the appliance away from any flame.
4. Keep flammable liquids away from the appliance.
5. Do not activate the ultrasound or the heating function if the tank does not contain cleaning liquid, to not damage the device.
6. Leave the filter at the bottom of the tank securely in place.

Maintenance

1. Periodically check the general condition of the device.
2. The device must be used by qualified personnel.
3. Check the condition of the electrical connections, the tank and the drain valve.
4. Do not allow dust to accumulate on the fan grille of the ultrasonic emitters.

Contents

- 1 Ultrasonic cleaner
- 1 Lid
- 1 Power supply cable
- 1 Basket
- 1 User manual

INSTRUCTIONS FOR USE

1. Check the condition of the plug and power cable for the ultrasonic bath.
2. Pour the cleaning solution into the stainless steel tank. Make sure the liquid level remains between 1/3 and 2/3 of the tank height.
3. Connect the device to the mains and turn on the switch (I).
4. For appliances equipped with the heating function, select the desired temperature by pressing the \uparrow or \downarrow arrows on the control panel. Press the «Heating» button to start the heating cycle (heating). A flashing dot on the temperature display indicates that the heating function is activated. The liquid temperature is displayed in real time and is regulated to achieve the desired temperature. After the selected temperature, the heating time may be long before reaching the set point.
5. Place the instruments to be cleaned in the basket, then place the basket in the tank resting on the edges.
6. Select the operating time by pressing the \uparrow or \downarrow arrows on the control panel (by pressing the average between 3 and 30 minutes).
7. Press the «ultrasound» button. The flashing dot on the minute display, ripples on the surface of the liquid and the sound of vibrations indicate that the cleaning cycle is in progress.
8. An audible signal indicates the end of the cleaning cycle.
9. Remove the basket from the tank and rinse the instruments.
10. Empty the tank manually by turning the drain valve.

Precautions for use

1. Read the instructions carefully before use.
2. Check the electrical connection before use.
3. Keep the appliance away from any flame.
4. Keep flammable liquids away from the appliance.
5. Do not activate the ultrasonic or heating function if the tank does not contain cleaning fluid to avoid Doing so will not damage the device.
6. Leave the filter at the bottom of the tank securely in place.

Maintenance

1. Periodically check the general condition of the device.
2. The device must be used by qualified personnel.
3. Check the condition of the electrical connections, the tank and the drain valve.
4. Do not allow dust to accumulate on the fan grille of the ultrasonic beams.

Contents

- 1 Ultrasonic cleaner
- 1 Low
- 1 Power supply cable
- 1 Basket
- 1 User manual

INSTRUCTIONS FOR USE

1. Check the condition of the plug and power cord of the ultrasonic bath.
2. Pour the cleaning solution into the stainless steel tank. Make sure the liquid level is between 1/3 and 2/3 remains the height of the tank.
3. Connect the appliance to the mains and turn on the switch (I).
4. For appliances equipped with the heating function, select the desired temperature by pressing the \uparrow or \downarrow arrows on the control panel. To start the heating cycle, press the «Heating» button (heating). A flashing point on the temperature display indicates that the heating function is enabled. The liquid temperature is displayed in real time and controlled to reach the desired temperature. After the temperature is selected, the heating time may take a long time to reach the set point.
5. Place the instruments to be cleaned in the basket and then place the basket in the tank, resting on the edges.
6. Select the operating time by pressing the \uparrow or \downarrow arrows on the control panel (by pressing the average to be pressed between 3 and 30 minutes).
7. Press the "ultrasound" button. The flashing dot on the minute display, the ripples on the liquid surface and the sound of vibration indicate that the cleaning cycle is in progress.
8. An audible signal indicates the end of the cleaning cycle.
9. Remove the basket from the tank and rinse the instruments.
10. Empty the tank manually by turning the drain valve.

Precautions for use

1. Read the instructions carefully before use.
2. Check the electrical connection before use.
3. Keep the appliance away from flames.
4. Keep flammable liquids away from the appliance.
5. Do not activate the ultrasonic or heating function when the tank does not contain cleaning fluid, otherwise it will damage the device.
6. Leave the filter in place at the bottom of the tank.

Maintenance

1. Regularly check the general condition of the appliance.
2. The appliance must be used by qualified personnel.
3. Check the condition of the electrical connections, the tank and the drain cock.
4. Make sure that no dust accumulates on the fan grille of the ultrasonic transmitters.

Contents

- 1 Ultrasonic cleaner
- 1 Lid
- 1 Power cord
- 1 Man
- 1 User Manual

INSTRUCTIONS FOR USE

1. Check that the plug and power cable of the ultrasonic bath are in good condition.
2. Pour the cleaning solution into the stainless steel tank. Make sure that the liquid level remains between 1/3 and 2/3 of the height of the tank.
3. Connect the appliance to the mains and turn the switch on (I).
4. For units equipped with a heating function, select the desired temperature by pressing the \uparrow or \downarrow arrows on the control panel. To start the heating cycle, press the "Heat" button.

(heating). A flashing dot on the temperature display indicates that the heating function is activated. The liquid temperature is displayed in real time and regulated to reach the desired temperature. Depending on the selected temperature, the heating time may be significant before reaching the set point.
5. Place the instruments to be cleaned in the basket then place the basket in the tank, resting on the edges.
6. Select the operating time by pressing the \uparrow or \downarrow arrows on the control panel (on average between 3 and 30 min).
7. Press the "ultrasonic" button. The flashing dot on the minute display, the ripples on the liquid surface and the sound of vibrations indicate that the cleaning cycle is in progress.
8. An audible signal indicates the end of the cleaning cycle.
9. Remove the basket from the tank and rinse the instruments.
10. Empty the tank manually by turning the drain valve.

Precautions for use

1. Read the instructions carefully before use.
2. Check the electrical connection before use.
3. Keep the appliance away from any flame.
4. Keep flammable liquids away from the appliance.
5. Do not turn on the ultrasound or the heating function if the tank does not contain any cleaning liquid, so as not to damage the device.
6. Leave the filter in place at the bottom of the tank.

Interview

1. Periodically check the general condition of the device.
2. The device must be used by qualified personnel.
3. Check the condition of the electrical connections, the tank and the drain valve.
4. Do not allow dust to accumulate on the fan grille of the ultrasonic transmitters.

Content

- 1 ultrasonic bath
- 1 cover
- 1 power cable
- 1 basket
- 1 user manual

INSTRUCTIONS FOR USE

- 1. Check the condition of the plug and power cord of the ultrasonic bath.
- 2. Pour the cleaning solution into the stainless steel tank. Ensure the liquid level remains between 1/3 and 2/3 of the tank's height.
- 3. Connect the device to the mains and turn on the switch (I).
- 4. For appliances equipped with a heating function, select the desired temperature by pressing the ÿ or ÿ arrows on the control panel. To start the heating cycle, press the "Heating" button.

(heating). A flashing dot on the temperature display indicates that the heating function is on. The liquid temperature is displayed in real time and regulated to reach the desired temperature. Once the temperature has been reached, the heating time may be long before reaching the set point.

- 5. Place the instruments to be cleaned in the basket and then place the basket in the tank, resting on the edges.
- 6. Select the operating time by pressing the ÿ or ÿ arrows on the control panel (pressing the average between 3 and 30min).
- 7. Press the "ultrasound" button. A flashing dot on the minute display, ripples on the liquid surface, and a vibrating sound indicate the cleaning cycle is in progress.
- 8. An audible signal indicates the end of the cleaning cycle.
- 9. Remove the basket from the tank and rinse the instruments.
- 10. Empty the tank manually by turning the drain valve.

Precautions for use

- 1. Read the instructions carefully before use.
- 2. Check the electrical connection before use.
- 3. Keep the appliance away from any flame.
- 4. Keep flammable liquids away from the appliance.
- 5. Do not activate the ultrasound or heating function if the tank does not contain cleaning fluid, to no do not damage the device.
- 6. Leave the filter at the bottom of the tank firmly in place.

Maintenance

- 1. Periodically check the overall condition of the device.
- 2. The device must be used by qualified personnel.
- 3. Check the condition of the electrical connections, the tank and the drain valve.
- 4. Do not allow dust to accumulate on the fan grill of the ultrasonic emitters.

Content

- 1 Ultrasonic cleaner
- 1 Way
- 1 Power cable
- 1 Road
- 1 User Manual

OPERATING INSTRUCTIONS

- 1. Check the condition of the plug and power cord of the ultrasonic bath.
- 2. Pour the cleaning solution into the stainless steel tank. Ensure the liquid level remains between 1/3 and 2/3 of the tank height.
- 3. Connect the device to the power supply and turn on the switch (I).
- 4. For devices with a heating function, select the desired temperature by pressing the ÿ or ÿ arrows on the control panel. To start the heating cycle, press the "Heating" button.
(Heating). A flashing dot in the temperature display indicates that the heating function is switched on.
The liquid temperature is displayed in real time and regulated to the desired temperature.
After the selected temperature, the heating time may take a long time before the setpoint is reached.
- 5. Place the instruments to be cleaned in the basket and then place the basket in the tank so that it rests on the edges.
- 6. Select the operating time by pressing the ÿ or ÿ arrow buttons on the control panel
(by pressing the average button between 3 and 30 minutes).
- 7. Press the "Ultrasonic" button. The flashing dot on the minute display, the ripples on the liquid surface, and the vibration sound indicate that the cleaning cycle is running.
- 8. An acoustic signal indicates the end of the cleaning cycle.
- 9. Remove the basket from the tank and rinse the instruments.
- 10. Empty the tank manually by turning the drain valve.

Precautions for use

- 1. Read the instructions carefully before use.
- 2. Check the electrical connection before use.
- 3. Keep the device away from any flame.
- 4. Keep flammable liquids away from the device.
- 5. Do not activate the ultrasound or heating function if the tank does not contain cleaning liquid to avoid damaging the device.
- 6. Leave the filter securely in place at the bottom of the tank.

maintenance

- 1. Check the general condition of the device regularly.
- 2. The device must be used by qualified personnel.
- 3. Check the condition of the electrical connections, tank and drain valve.
- 4. Make sure that no dust accumulates on the fan grille of the ultrasonic transmitters.

Contents

- 1 ultrasonic cleaner
- 1 lid
- 1 power supply cable
- 1 basket
- 1 user manual